

Laser Compatibility Chart for Optical Integrity ScopeSafe™ Holmium Laser Fibers

Rev B effective 22 Mar 2017

Brand Name	Manufacturer	Power (Watts)	Connector Type
Versa Pulse	Lumenis/Coherent	100 W	Standard SMA 905
Versa Pulse	Lumenis/Coherent	60 W	Standard SMA 905
Versa Pulse	Lumenis/Coherent	20 W	Standard SMA 905
DUR-HL20	Dornier/ACMI	20 W	Standard SMA 905
StoneLight	AMS	15 W	Standard SMA905 +metal adapter
i) MBL 70/20 RevoLix DUO ii) MCL 70W RevoLix 70 iii) MCL 120W RevoLix 120 iv) LCL 15W RevoLix Jr v) MPL 80W Sphinx 80 vi) MPL 45WLitho Sphinx Litho	Lisa Laser	See brand	Standard SMA 905 with slim sleeve only
Themis	Direx	30W	Standard SMA 905
Medilas H20	Wavelight	20W	Standard SMA 905
Quanta LithO	Quanta	30W	Standard SMA 905
Odyssey	Convergent	30W	Special fiber required
CALCULASE I	KARL STORZ	10 W	Standard SMA 905 with slim sleeve only
Omni Pulse MAX	Trimedyne	80 W	Special fiber required

Note: This list is subject to change as new lasers and new fiber models are introduced. Also, not all fiber models are compatible with all lasers, so if a laser is not on this list, contact us and we will work with you to ensure that you get a compatible fiber.