



Compatible for Use with Aura Laser Systems*

Product Features:

- Fiber attaches directly to laser without the need for adaptor
- Fiber tip is pre-stripped for procedural efficiency
- LUER connection secures fiber to hoop dispenser

Cross Reference Information:

Core Size	BSC Part #	Optical Integrity Part #
0.4mm/400µm	10-0622	11638
0.6mm/600µm	10-0612	11639

Order Information:

Part #	Description
11638	400µm Aura-compatible laser fiber, single use, sterile
11639	600µm Aura-compatible laser fiber, single use, sterile